



LIQUID URINE CONTROL - LEVEL 2 (LIQ URN CONTROL 2)

Cat. No. UC5074 / Lot No. 1237UC Size: 10 x 10 ml / Expiry: 2024/08/28

| Analyte / Instrument Group | Unit | Target | Range | | 1 - SD | Method |
|-----------------------------------|-------|--------|--------|---------|---------|---------------------------------|
| | | | Low | High | | |
| Creatinine – Mindray BS Series | mg/dL | 66,17 | 51,61 | 80,73 | 7,279 | Creatinine PAP |
| Glucose – Mindray BS Series | mg/dL | 22,11 | 17,24 | 26,97 | 2,432 | Glucose Oxidase |
| Microalbumin – Mindray BS Series | mg/L | 25,19 | 19,65 | 30,73 | 2,771 | Immunturbidimetric |
| Phosphate – Mindray BS Series | mg/dL | 23,34 | 18,20 | 28,47 | 2,517 | Phosphomolibdate UV |
| Total Protein – Mindray BS Series | mg/dL | 13,89 | 10,83 | 16,95 | 1,528 | Pyrogallol Red |
| Urea – Mindray BS Series | mg/dL | 918,46 | 716,40 | 1120,52 | 101,031 | Urease Kinetic |
| Uric Acid – Mindray BS Series | mg/dL | 8,99 | 7,02 | 10,97 | 0,989 | Uricate Peroxidase no Ascorbate |

| | | Low | High | 1 - SD | Method |
|-------|-------|------------|-------------|---------------|---------------|
| ug/dL | 21,14 | | | | |
| pg/mL | | | | | |
| ng/mL | | | | | |
| ng/mL | | | | | |

| | | | | | |
|-----------------|-------|--------|-------|-------|--|
| | pg/mL | 559,20 | | 3,146 | |
| 25-OH Vitamin D | ng/mL | 28,60 | 22,31 | | |