



LIQUID URINE CONTROL - LEVEL 3 (LIQ URN CONTROL 3)

Cat. No. UC5075 / Lot No. 1238UC Size: 10 x 10 ml / Expiry: 2024/08/28

Analyte / Instrument Group	Unit	Target	Range		1 - SD	Method
			Low	High		
Creatinine – Mindray BS Series	mg/dL	165,17	128,84	201,51	18,169	Creatinine PAP
Glucose – Mindray BS Series	mg/dL	269,34	210,09	328,59	29,627	Glucose Oxidase
Microalbumin – Mindray BS Series	mg/L	59,01	46,03	71,99	6,491	Immunoturbidimetric
Phosphate – Mindray BS Series	mg/dL	37,49	29,24	45,74	4,124	Phosphomolibdate UV
Total Protein – Mindray BS Series	mg/dL	68,54	53,46	83,62	7,540	Pyrogallol Red
Urea – Mindray BS Series	mg/dL	1728,05	1347,88	2108,23	190,086	Urease Kinetic
Uric Acid – Mindray BS Series	mg/dL	13,85	10,81	16,90	1,524	Uricate Peroxidase no Ascorbate

